## Instructions for use

en

### PRONOLIS HD

Intra-articular sodium hyaluronate

#### 1. Description

The prefilled syringe contains a hydrogel made from highly purified sodium hyaluronate obtained by bacterial fermentation. This gel is used to improve viscosity of intra-articular connective tissues. It has been subjected to moist heat sterilization for injectable products and is non-pyrogenic.

Depending on the product version and concentration, 1 ml aqueous, physiological phosphate-buffered solution contains the following quantity of sodium hyaluronate: 10 mg, 16 mg, 22 mg or 25 mg.

The sodium hyaluronate hydrogel is therefore contained in a weakly saline, sterile, non-pyrogenic isotonic phosphate-buffered saline at pH 7.2 (6.8-7.4) in a volume sufficient for injection. Depending on the product version, full volume is 1.4-8 ml.

Chemical composition (w/v) as mg of substance per 1 ml hydrogel: see table.

#### 2. Properties and Mode of Action

A substantial component of healthy synovial fluid is hyaluronic acid, a ubiquitous glycosaminoglycan. This natural biopolymer maintains viscoelasticity, and its lubricating and shock-absorbing properties enable pain-free physiological joint movements. It also participates in cartilage nutrition. The sodium hyaluronate contained in the product is a salt of hyaluronic acid.

#### 3. Indications

Pain and decreased articular mobility associated with degenerative lesions of the knee and other synovial joints (hip, ankle, shoulder, elbow, hand), including osteoarthritis.

#### 4. Contraindications

The product should not be used in patients

- with known hypersensitivity to any of the components of the product.

- with joint inflammation.

- with septic arthritis.

- with skin infections or diseases at the injection site.

- with natural or drug-induced clotting deficiency, for example due to hemophilia or use of anticoagulants such as Marcumar (phenprocoumon) or Coumadin (warfarin).

Children, pregnant or nursing women should not be treated with hyaluronic acid as there are no clinical data available on its use in these patients.

#### 5. Precautions

Aspects of arthritis is a potential serious side effect of treatment, please observe all the precautionary measures usually required for surgical procedures. The hydrogel should be injected into the joint cavity. Intra-vascular and intra-tissue injection should be avoided. It must not be resterilized, resealed or reused. Do not sterilize the solution, as this can cause changes to the product characteristics. Do not use the product if there is any damage to the ready-to-use-syringe or sterile packaging. The product is to be used for one patient only during a single appointment.

#### 6. Potential Side Effects

Biocompatibility testing of the product did not show undesirable effects. Transient local reactions, such as redness, swelling, sensation of heat and pain, may occur after intra-articular injection. These effects can be reduced by cooling the area with ice for 5 to 10 minutes after the injection. During the first days of treatment, concurrent oral intake of analgesics/anti-inflammatory drugs (NSAIDs) may prove beneficial. Treatment must be stopped if the patient experiences intense pain or inflammation. Adverse reactions will be treated under the supervision of the treating physician. Strenuous physical activity should be avoided directly after the injection. All instructions provided by the physician must be followed.

The healthcare professional should advise the patient

- of any potential adverse events related to the product.

- that the patient must report any adverse event or complication to a physician.

Serious events that have occurred in connection with the use of the product must be reported to the manufacturer and to the competent national authority in which the physician and/or patient is established.

#### 7. Drug and Chemical Interactions

To date, there are no data about incompatibility with other intra-articular injections. There exists a chemical incompatibility between sodium hyaluronate and quaternary ammonium compounds such as benzalkonium chloride that may lead to crystal precipitation. Therefore, this product should never come in contact with surgical instruments or other products that contain them or were rinsed with solutions containing quaternary ammonium compounds as biocidal preservatives.

#### 8. Dosage, Mode and Duration of Treatment

To reduce risk of inflammation the following steps should be followed. Take the prefilled syringe from the packet. The protective closure of the Luer lock is removed. A suitable sterile needle is screwed onto the syringe. Use of a single-neededle is recommended. Remove air from the syringe before injection. Direct the designated area of the skin with a suitable topical antiseptic so that the product can be administered under the best conditions.

If there are signs of intra-articular swelling, fluid should be aspirated and tested for a bacterial etiology before the product is injected. The hydrogel is injected intra-articularly 1 to 5 times at intervals determined by the physician. The beneficial effects of the treatment last at least 6 months. Treatment cycles may be repeated, if required. Depending on the full volume, a patient may have several joints treated at once. Discard any unused product.

#### 9. Packaging

Box containing a sterile prefilled syringe with 1-4.8 ml.

#### 10. Shelf Life

Do not use the product after the expiry date. The expiry date (year/month) is printed on the syringe blister pack and box.

#### 11. Storage Conditions

The product must be stored dry between 2 °C and 25 °C in the original packaging and protected from light and impact damage. Avoid freezing or impact. Please observe the symbols on the package. Store this product in a location inaccessible to children.

#### 12. User Information

If a product from damaged packaging is used, or a product reused on other patients, or on the same patient at a later time, the intended properties and sterility of the product cannot be guaranteed. Potential material contamination and transfer of pathogens from the general environment, the clinical environment or (other) patients could lead to health complications. Resealing or resterilization is not a permanent or guaranteed method of restoring the product to a useable condition.

#### 13. Recommended user type

For exclusive use by physicians.

#### 14. Disposal

The syringe content is non-toxic and non-flammable. Unused syringes and their contents are not infectious and may be disposed in compliance with national and local regulations, before or after their expiry date. Used syringes and needles must be disposed of as epidemiologically hazardous waste according to national and local regulations for their safe handling and disposal.

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