Natriumhyaluronat intraartikulär

Formel	PRONOLIS	PRONOLIS ^{HD} mini	PRONOLIS ^{HD} one	PRONOLIS ^{HD} mono
1 ml enthält:	1,0 %	1,6 %	2,2 %	2,5 %
Natriumhyaluronat	10,000 mg	16,000 mg	22,000 mg	25,000 mg
Natriumchlorid	8,500 mg	8,500 mg	8,500 mg	6,500 mg
Dinatriumhydrogenphosphat 2H ₂ O	0,563 mg	0,563 mg	0,563 mg	0,563 mg
Natriumdihydrogenphosphat 2H ₂ O	0,045 mg	0,045 mg	0,045 mg	0,045 mg
Wasser für Injektionszwecke	q.s.	q.s.	q.s.	q.s.

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Origin
Sodium hyaluronate is obtained by bacterial fermentation and is used as a viscoelastic gel to supplement the intra-articular connective tissue. It is sterile and apyrogen.

Composition
1 ml of solution contains 10 mg, 16 mg, 22 mg, 25 mg of sodium hyaluronate with fermentative origin, a low phosphate tampon sterile, isotonic and apyrogen, pH 7.2 and water for injection.

Indications
It is injected into the knee to relieve pain and decreased articular mobility associated with degenerative knee lesions and other synovial joints.

Contra-Indications
The product should not be used for patients
- with a known hypersensitivity reaction to any of the components of the product
- with septic arthritis
- with infections of the skin or dermatological diseases at the injection side
- given coagulation inhibitors such as Marcumar

Cautions
Since septic arthritis is a serious side effect, please observe all usual precautions for surgical interventions.
The product has to be used before the expiry date indicated on the packaging.
No sign of bioincompatibility has been detected throughout biocompatibility testing of the product. Some local transient reactions may occur after intra-articular injection, such as redness, swelling, heat and pain. The local application of ice for 5-10 minutes may reduce these effects.
No known incompatibility with other intra-articular injectable products has been reported up to now.
Simultaneous intake of analgesics or anti inflammatory drugs by oral route may be useful during the first days after treatment.
The product must not be used in association with quaternary ammonium compounds.
The product must not be resterilised or reuse because this could lead to a risk of infection and/or cross-infection and/or malfunction of the product.

Administration and dosage
The syringe is prefilled, the Luer lock closure is removed and a suitable sharp edged tip sterile needle is screwed on the Luer lock. Use of a single use needle is recommended. On desire, needles are included. Remove air bubbles from the syringe before injection.
The skin is disinfected with local antiseptic.
The product is injected in the joint, once a week, 1-5 times. The beneficial effects of 1-5 sequential treatments are maintained during at least 6 months. In the case of intra articular swelling, with fluid should be removed by aspiration before the hyaluronic acid injection.
Depending on the size of the joint: 1, 3 or 5 injections of 1 ml or 2 ml may be administered intra-articular.
The prefilled syringe of 1 ml is the volume adapted for small joints.
The prefilled syringe of 2 ml is the volume adapted for the knee, hip, shoulder and ankle.
The prefilled syringe of 4.8 ml is the volume adapted for the knee only.
The recommended treatment regimen with the mono dose of 4.8 ml for the knee is one injection, which may be repeated 6 months later based on the clinical experience and the pain relief obtained by the patient.

Properties and mode of action
The synovial fluid contains physiologically hyaluronic acid which is responsible for viscoelastic properties allowing, through lubricant properties, painless joint movements. It also participates to the nutrition of the cartilage.
Intra-articular supplementation of the synovial fluid by hyaluronic acid injections improves the viscoelastic properties of the synovial fluid. After 1-5 intra-articular injections, with joint mobility is improved, and the pain associated with degenerative lesions is relieved.

Packaging
a) 1 box containing a prefilled syringe à 1 ml or 2 ml.
b) 1 box containing a prefilled syringe of 4.8 ml containing 120 mg sodium hyaluronate within sterile packaging.

Type of user recommended
Exclusively medical doctors

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ALBOMED® GmbH · Hildebrandstrasse 11 · 90592 Schwarzenbruck, Germany